

Xu lab protocol and instructions

Paraformaldehyde (PFA) solutions

We make 16% PFA stock solution, and make 4% PFA for animal perfusion and slice fixation.

Materials Needed:

-Weigh boats

-Thermometer

-Filter paper

-Funnel

-2 beakers w/H₂O

-Big scooper

-Graduated cylinder

-pH indicator strips

-Scale

-- hot and stirring plate

-- 1 liter volumetric flask

-KOH, NaOH for getting to final **pH 7.2**

We will make PFA solution in the mouse vivarium (Room 338) procedure room with a chemical hood.

16% PFA

- **Measure 160 g PFA for making one liter of 16% PFA**
- Heat ~700-800mL water to no more than 60°C – Make sure not to boil!

Xu lab protocol and instructions

- Add PFA and stir with stirring
- Add 2-3 drops of 10N NaOH
- Allow sufficient time for PFA to dissolve
 - If takes a long time, add more NaOH
- Filter solution using a fine or medium porosity filter
- Adjust to pH 7.2 adding drops of HCl or NaOH, depending on if pH needs to be increased (NaOH) or decreased (HCl). // using pH indicator strips
- Top off with water to 1L in a 1 liter volumetric flask

4% PFA

Measure each solution below with an appropriate graduated cylinder

Mix thoroughly with the appropriate proportions:

Solution A: 1 part 50mL 5mL 2mL

Solution B: 4 parts 200mL 20mL 8mL

16% PFA: 2.5 parts 125mL 12.5mL 5mL

ddH₂O: 2.5 parts 125mL 12.5mL 5mL

500mL 50mL 20mL